

The revolution in Dry Eye analysis.

The new solution for smart Dry Eye analysis.

Dear Eye Care specialist,

are noticing a change in their people's quality of life. visual comfort. The habits of modern life, with the intensive The concept of tearcheck® can feel dry and tired.

most common eye disease worldwide and has a strong impact on people's lives and visual comfort.

The earlier Dry Eye cases are evaluated, the better Dry Eyes can be managed. The easier Dry Eye analysis can be done, the more frequently patients' eyes can be checked. The On behalf of the entire R&D team, the changes to their visual innovation come true

Taking care of vision comfort comfort, the easier it will be is an exciting mission, knowing to provide working and lasting that more and more people solutions to finally improve

use of digital devices is only was driven by the idea of one of the reasons why eyes creating a new standard for smart, reliable and efficient The condition of Dry Eyes Dry Eye analysis, which allows is affecting all generations, to evaluate patients' visual with a rapid increase in the comfort fast, easily and fully younger generation. Dry Eve automatically, tearcheck® has Disease became the second been designed based on the latest technology.

> Be inspired by our innovation and discover the difference in daily Dry Eye Management.

Yours, Petra BRUCKMUELLER CEO ESW vision

simpler it is to show patients who made this break-through



Visionary technology

tearcheck[®] allows **quick and** easy examinations, as in total 9 exams can be done, which provide **break-through Dry** Eye analysis in less than 10 minutes, report included. Two high-resolution cameras inside provide perfect capture Also patented is the exam used to perform the different exams within a self-explaining process, guiding you through.

A patented new way to evaluate tearfilm stability on the whole cornea, by using a line mask, is provided by two exams, which are done within one recording (TFSE[®]/NIBUT).

of images and videos. White, to evaluate the degree of blue and infrared light are inflammation on the ocular surface (OSIE®).





Evaluation

OSDI Ocular Surface Disease Index



Tear Meniscus

for smart analysis




OSIE® Ocular Surface Inflammatory Evaluation



Dynamic **Breakup** Time









Demodex Lash base alteration




Eye Redness



Abortive Blinking Eyelid Tonicity

break-through Dry Eye analysis in less than 10 minutes



A revolution is born.

Taking care of Dry Eye tearcheck[®] supports multipatients and their visual user management. Every user **comfort requires educating** has their own profile, with their them about risks, causes, own patient database. signs and impact.

different climates, air pollution, it on the table, connect the frequent use of digital screens tablet and start. and aging are well-known reasons for more and more tearcheck® provides fully people claiming to have less automatic analysis for timevisual comfort.

tearcheck[®] is a stand-alone device with an integrated display. This screen allows the tearcheck® can be used user to view all the acquisitions wherever it fits your workand examinations in real time. **place.** No dark room required, The user can also wirelessly no computer connection of his choice (tablet, smart- connection needed. phone, PC, etc.) in order to duplicate the tearcheck[®] screen. And last but not least:

patient management system which can be shown directly inside: every exam can be on the screen by selecting stored and shown to the concerned dates and types of patient at any time to follow exams, exported via USB in pdfup the progress or changes in format, or printed wirelessly. patients' eyes.

tearcheck[®] is delivered com-The modern way of life, pletely installed. Just mount

> efficient ease-of-use. It simply makes it fun to get things started.

connect an additional device required, no external cable


tearcheck® results are tearcheck[®] comes with a presented in a final report As you prefer ...

Be inspired by the innovation



... to create a disruptive product fully dedicated to Dry Eye syndrome and capable of constantly evolving, that was my intention behind developing tearcheck[®].

> Dr. Yves Vincent **BROTTIER** Creator of tearcheck®



TFSE® **Tear Film Stability Evaluation**

Makes it possible to assess the micro-deformations appearing on the surface of the tear film. These micro-deformations reflect instability of the tear film. They are presented in number and in intensity (score and diagram). And recorded within 10 seconds.

TFSE[®]/ NIBUT

A new kind of exam to evaluate the surface of the tear film three times per second. Specifically designed optimised horizontal lines replace the traditional circular pattern.




tearcheck[®] line mask



NIBUT

Non-Invasive **Breakup-Time**

Shows the **ruptures in the tear film** already from two seconds on, based on highly sensitive detection. Recorded within TFSE[®] exam.







Meibography 2D imaging / 3D rendering

Makes it possible to visualize the Meibomian glands in 2D imaging and 3D rendering. The results help to evaluate the rate of gland loss (in %) and the **morphology** of the glands present.

A special eyelid flipping tool is used to allow touch-free examination.



tearcheck[®] eyelid flipping tool



Eye Redness

This examination captures images of the bulbar and palpebral conjunctivae. It makes it possible to **assess** hyperemia induced by the inflammatory phenomena consecutive to Dry Eye syndrome and to monitor its progress over the treatments delivered to the patient.

Demodex

The Demodex exam consists of an enlarged image capture of the base of the eyelashes. The magnification factor of the image makes it possible to find and visualize the signs of the presence of the parasite.







Abortive Blinking

Within 60 seconds, all blinks are recorded for automatic **detection** of incomplete blinking.





OSIE®

Ocular Surface Inflammatory Evaluation

The principle behind this examination resides in the fact that fluorescein will remain "attached" to areas of the eyeball exhibiting surface alteration in connection with an **inflammatory process** in the case of Dry Eye syndrome.





D-BUT

Dynamic Breakup Time

Additional exam to measure **tear film breakup time** with a drop of fluorescein.





OSDI Ocular Surface

Ocular Surface Disease Index

A questionnaire intended to be filled out by the patient, reflecting their subjective feeling and opinion on their **visual comfort** in their daily life.



Is evaluated in the lower eye lid and is showing, if there is sufficient **tear production** before the next eye blinking.

C OSDI									
General well-being:	Yes			No					
Do you feel good today?		\checkmark							
Are you under medical treatment with an ophtalmologist?				\checkmark					
Are you especially sensitives to light?				\checkmark					
Are you pregnant?	\checkmark								
Analysis of eye fitness:				2					
Did you experience the following phenomena last week?	Never	Sometimes	Regulary		Often	Always			
Light-sensitive eyes			\checkmark						
Sandy feeling									
Sensitives and irritated eyes									
Blurred vision									

13

te	eck	Report History								
Ĺ	-	Date	Exam	Ô	P					
◎ ←		2000-03-03 01:29	NIBUT			\checkmark				
		2000-03-03 00:46	NIBUT							
		2000-03-03 01:29	TFSE			\checkmark				
		2000-03-03 00:46	TFSE							
		2000-03-03 01:17	OSIE							
_	_	2000-03-03 01:15	OSIE							
		2000-03-03 01:34	Demodex							
25	32 22	2000-03-03 01:35	Eye Redness							







USB export

wireless print view on the screen



Reports

All reports per patient will be stored. This makes it possible to compare different results from same examination or regard all selected exams. Report history can be filtered by date or exam.

Discover the difference

tearcheck® was designed to tearcheck® software is set up create a new standard for to allow ongoing developsmart, efficient and reliable ment, providing innovative Dry Eye analysis. The ease- opportunities of use during its of-use combined with highly whole lifetime period. innovative technology inside to provide accurate, brilliant Smart connectivity will open and representative images, new possibilities of interaction will bring a new comfort of between Dry Eye specialists examination in daily Dry Eye and patients. Management.


Ease of Use

Discover the comfortable, selfoperation exaplantory of tearcheck[®]. Review the results directly with the patient on the wirelessly connected tearcheck® tablet.



15






Marketing support

exciting marketing support direct from supplier for your practice. tearcheck® comes with a marketing package your user experience. available on demand to make the best out of your Dry Eye Management education and advertising. Marketing support in the form of manuals, bro-

With tearcheck[®], benefit from chures, POS-collateral, social media posts, press and online media are there to add value to your practice and enhance

> Visit www.esw-vision.com, or contact our marketing team (marketing@esw-vision.com) for more information!





tearcheck® smart & complete

Technical Data

Camera	2 high-resolution cameras			
Storage space	200 GB			
Light sources	white light / blue light (465nm) / near infrared light (850nm)			
Certification	CE class I / ISO 13485			
Conformity	93/42/EEC, directive Annex II, article 3			
Light emission	no risk group, according IEC60601-2-57			
RF emission	2.5 GHz			
Contra-indications	-			
Shipping dimensions (LxWxH)	815×520×355 mm			
Shipping weight	12.5 kg			
Energy Source	100-240V / 50-60Hz external power supply			
Dimensions (LxWxH)	240×195×190mm (active part only)			

Designed for innovative Dry Eye analysis

Easy I smart I unique

Discover the complete solution for analysis and treatment.

Dry Eye Management. One step ahead.

www.esw-vision.com



Tune in for more insights.

www.esw-vision.com

tearcheck®, OSIE® and TFSE® are registered brand marks of E-Swin, medical device, class I Ref. Nr. M040GB003B____ rev. B 2021/04/19