## optimis

Cataract and Glaucoma Laser







ootimis

glaute optimis!

Quantel Medical's cutting edge laser technology has brought constant and multiple innovations to ophthalmologists worldwide since **1993**.

Optimis II is a 1064nm Nd:YAG laser integrated into a premium diagnostic slit lamp.

## Optimis II from Quantel Medical.

Advanced YAG Laser Technology. Optimis II integrates a YAG laser with diagnostic slit lamp offering a compact, all-in-one design. The Optimis II combines premium optics and a powerful laser cavity to provide optimum photodisruption for capsulotomy and iridotomy procedures.

### **Diagnostic Slit Lamp**

- Convergent Optical Path
- Optimized for viewing the anterior and posterior segment

## Advanced YAG Technology

- Highly efficient, powerful and reliable
- Precise two-point aiming system
- Ajustable posterior offset

### **Advanced Design**

- Integrated platform
- Compact footprint
- Combines with a green laser


## **YAG Mode**

The YAG mode delivers high-performance photodisruption for capsulotomy and peripheral iridotomy. The gaussian laser beam profile allows for precise delivery, only requiring minimum energy levels, which avoids side effects such as lens pitting.

Precise tissue targeting is achieved with the advanced two-point aiming beam, allowing for fast and accurate targeting of the capsule; and the adjustable offset feature, allowing for easy focus adjustment behind or at the point of disruption.

### Precise and Powerful





# optimis



## TECHNICAL SPECIFICATIONS

#### MH-VAG

Laser source: Q-Switched Nd:YAG
Wavelength: 1064nm

Energy levels: 0.5mJ to 10mJ (1 pulse)

Pulse duration: 4n

Spot Size: 10 µm (86.4% energy in diameter)

8 μm FWHM Repetition rate: Up to 1.5Hz

Offset: +30µm, +100µm, +200µm
Pulse mode: 1, 2 and 3 pulses per shot
Aiming beam: Red 650 nm, adjustable intensity

CSO 980

Magnification: 3 positions: 10x, 16x, 25x

Air cooled

Size (laser only): 67 (H) x 41 (W) x 50 (D) cm 24.4" (H) x 16.1" (W) x 19.7" (D)

Weight (laser only): 24.6kg 54.3lbs

Power requirements: 100 – 240VAC, 50/60Hz

### **AVAILABLE MODELS**

Slit Lamp Type:

Cooling:

Optimis II

Optimis II and Vitra laser platform

Optimis II and Vitra Multispot laser platform

Specifications are subject to change without notice.

©2014.Quantel Medical, Optimis II is a registered trademark of Quantel Medical. All rights reserved.

www.quantel-medical.com

Designed and Manufactured by Quantel Medical in France







### Headquarters

Quantel Medical
11, rue du Bois Joli - CS40015
63808 Cournon d'Auvergne – FRANCE
Tel: +33 (0)4 73 745 745

Fax: +33 (0)4 73 745 700 E-mail: contact@quantel-medical.fr

### North America

Quantel USA 601 Haggerty Lane Bozeman, MT 59715 – USA Tel: +1 877 782 6835

Fax: +1 406 522 2005

E-mail: info@quantelmedical.com

**Representative Offices** Thailand, Chiang Mai Brazil, Rio De Janeiro

